### BHL-Multiflow®

**Contrast Media Administration System in T** - REF: 4000103-B  
**Contrast Media Administration System in T** - REF: 4001111  
**Contrast Media Administration System in T** - REF: 4001111C

#### Indications


BHL- Multiflow® are destined for Radiology procedures, CT and MRI, with administration of Contrast Media. The devices identified are suitable for 12h of use or 50 injections of contrast media, whatever is achieved first, when connected to the *SingleFlow®* patient lines (BHL-SVM, BHL-DVM, BHL-SVC and BHL-DVC) and under aseptic conditions. Asepsis procedures must be strictly complied to reduce the risk of contamination. The *Contrast Media Administration System in T* must be used only by trained personnel to administer contrast media to a patient. The *Contrast Media Administration System in T* is compatible with all the syringes equipped with luer connections. The *Contrast Media Administration System in T* connected to a patient line *Singleflow®* allows repeated injections of contrast media. The patient line *Singleflow®* must be changed for each patient. Medical personnel are responsible for preventing any cross contamination. The user is responsible for ensuring that if any other devices are used with the *Contrast Media Administration System in T*, they are compatible and adapted for repeated use of the *Contrast Media Administration System in T*. Information supplied with the medicines to be used in *Contrast Media Administration System in T* must be taken in consideration before any procedure. The products are Luer connections and completely phthalate free. The safety use of the *Contrast Media Administration System in T* is assured if the protocol described is followed.

#### Protocol

1. Connect the spike of the drip chamber with blue stripe tube to the bottle/bag with saline solution. Only for REF: 4000103B.
  2. Connect the spike or the drip chamber spike with non stripe tube to the bottle/bag with contrast media.
  3. Connect the T connectors on the syringes of the injector.
  4. Connect the patient line *Singleflow®* to the *Contrast Media Administration System in T*.
  5. Purge the system and ensure the air is removed from the tubing system.
  6. Assure that all the connections are correctly executed and start the radiologic procedure.
  7. When the procedure is finished disconnect the used *Singleflow®* patient line from the *Contrast Media Administration System in T*.
  8. Connect a new *Singleflow®* patient line into the *Contrast Media Administration System in T*. Leave the cap in the terminal part of the patient line.
- Only remove the cap when it's to be connected to the catheter to be introduced in the patient.

#### Warnings!

- If residues are observed at any point on the tubing system, including the particle filter, the *Contrast Media Administration System in T* must be replaced.
- Do not exceed the pressure of 300Psi. The use of a greater pressure or occlusions in the fluid path might result into leaks or ruptures.
- Needles or catheters with reduced flow influence the working pressures of the system.
- It's recommended to use multipatient syringes with the *Contrast Media Administration System in T* system.
- Before start a new produce make sure the patient line is changed.
- When it's necessary to change the automatic injector syringes the *Contrast Media Administration System in T* must also be changed.
- After use devices must be treated as hospital waste.
- The manufacturer is not responsible for damages that may arise from the device misuse or for the non accomplishment of the Instructions for use.
- Defined storage conditions must be kept and assured.

 Manufacturer	 Sterilized using ethylene oxide	 Liquid filter with pore size
 Use by date	 Do not use if package is damaged and consult instructions for use	 One-way valve
 Batch code	 Single sterile barrier system with protective packaging outside	 Consult instructions for use or consult electronic instructions for use
 Catalogue number	 Keep dry	 Non-Pyrogenic
 Distributor	 Keep away from sunlight	 Medical device
 Country of manufacturer	 Caution	
 Do not resterilize		

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




## BHL-Multiflow®

**Contrast Media Administration System in T - REF: 4000103-B**

**Contrast Media Administration System in T - REF: 4001111**

**Contrast Media Administration System in T - REF: 4001111C**

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