BHL-Hiflow®

High Pressure lines with male/female connector and valve; High Pressure lines with male/male connector and valve; High Pressure lines with male/female connector; High Pressure lines with male/male connector; High Pressure lines with female/male rotative connector; High Pressure lines with male/male rotative connector; High Pressure lines with male/male connector and stopcock; High Pressure lines with male/female connector and stopcock; High Pressure lines with male/male connector, stopcock and protective mango;

Indications

High Pressure Lines and Accessories are intended to be used for administration of contrast media into the patient body in medical procedures that require de administration of fluids under High Pressure. They are intended to be connected to the syringes of an automatic injector to administrate the fluid, contrast media, under pressure into the patient body. One side of the device connects to the syringe and the other end connects to the catheter placed on the patient. The fluid is then injected from the syringe, through the High Pressure line into the catheter. The luer connections of the device make them compatible with the syringes of the automatic injector and with other medical devices with luer connections. The High Pressure Lines must be used only by trained personnel. The High Pressure Lines are compatible with all the syringes equipped with luer connections. The user is responsible for ensuring that if any other devices are used with the High Pressure Lines, they are compatible. The High Pressure Lines are Luer connections and completely phthalate and latex free. The safety use of the devices is assured if the protocol described is followed.

Protocol

1. Connect one end of the device to the automatic injector syringe;
2. Purge the High Pressure line and ensure the air is removed from the tubing system.
3. Connect the other end of the device to the catheter placed on the patient. Only remove the cap when the High pressure line it's to be connected to the catheter to be introduced into the patient.
4. Start the procedure.
5. When the procedure is finished disconnect the High pressure line from the syringe of the automatic injector and from the catheter placed on the patient.

Warnings!

- Devices are suitable to use under pressures up to 1200 Psi. The use of a greater pressure or occlusions in the fluid path might result into leaks or ruptures.
- Needles or catheters with reduced flow influence the working pressures of the devices
- All the Medical devices part of this IFU are provided sterile, single-use condition, and are intended to be used only once for a single patient.
- Side effects described for this type of product are cross-contamination, which can result from multiple use of the product, reduction of professional hygiene and environmental conditions, and use of the product with a damaged packaging. Assepsis procedures must be strictly complied to reduce the risk of contamination.
- After use devices must be treated as hospital waste.
- The manufacturer is not responsible for damages that may arise from the device misuse or for the non accomplishment of the Instructions for use.
- Defined storage conditions must be kept and assured.

Manufacturer	Sterilized using ethylene oxide	Caution
Use by date	Do not use if package is damaged and consult instructions for use	Do not re-use
Batch code	Single sterile barrier system with protective packaging outside	Consult instructions for use or consult electronic instructions for use
Catalogue number	Keep dry	Non-Pyrogenic
Distributor	Keep away from sunlight	MD Medical device
Country of manufacturer		
Do not resterilize		

BHL-Hiflow®

High Pressure lines with male/female connector and valve; High Pressure lines with male/male connector and valve; High Pressure lines with male/female connector; High Pressure lines with male/male connector; High Pressure lines with female/male rotative connector; High Pressure lines with male/male rotative connector; High Pressure lines with male/male connector and stopcock; High Pressure lines with male/female connector and stopcock; High Pressure lines with male/male connector, stopcock and protective mango;

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